



Test kit color disc chlorine dioxide, model CLD-2

Številka izdelka:

2842800

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan

Nearno



Monitor ClO₂ levels on-site with the NEW Chlorine Dioxide Test Kit. The Chlorine Dioxide Test Kit uses the reliable and EPA-accepted DPD method to monitor ClO₂ concentration, while eliminating chlorine interference through the addition of glycine. With 100 tests per kit and available range of 0-2 mg/L and resolution of 0.1 mg/L, operators can accurately monitor around the 0.8 mg/L regulated limit.

Simple, cost-effective method for maintaining optimum chlorine dioxide levels in your plant

USEPA-accepted for waste water and potable water

Approved by some states for potable water under specific conditions (check with your regulatory agency)

Specifikacije

Ime metode:	Color Disc/DPD
Merilno območje:	0 - 2.0 ppm ClO ₂
Model:	CLD-2
Parameter:	Chlorine Dioxide
Platforma:	Color Disc
Smallest increment:	0.1
Standardno merilno območje:	0 - 2.0 ppm ClO ₂
Število testov:	100
Test Kit model:	CLD-2 ²
Vrsta primera:	D